



**Totaalleverancier
van lasers en medische
technologieën**



Laservision Instruments BV
+31 (0) 529 428 000

www.laservision.nl
info@laservision.nl

INDIBA®
revitalizing lives



ONA



- ① Hydrodermabrasion
Head Cleansing
(ONA handpiece)
► SOLUTIONS
- ② Electroporation + LED
Head Infuse
(ONA EP handpiece)
► ONA MOISTURIZING SERUM
- ③ Low-frequency electrotherapy + Massage
Head Workout
(ONA Y handpiece)
► ONA GEL

ONA breaks the rules.

The INDIBA® facial device combines hydrodermabrasion (double peel), electroporation, LED, and low-frequency electrotherapy to elevate your skin to a new level of health and purity.

ONA is pleasure, sensations, purity, and perfection. Optimize and enhance your skin with the ONA trend.

Key Benefits

- Exfoliates, cleans, and removes impurities
- Deeply moisturizes
- Brightens and unifies skin tone
- Increases muscle tone, helping to improve facial contours
- Improves the overall appearance of the skin

Technology

At INDIBA®, innovating with different technologies is key to satisfying the needs of our clients.

ONA is designed to work with different technologies:

○ CHEMICAL PEELING

Eliminates sebum and skin impurities, all while providing anti-aging active ingredients through our ON solutions.

○ MECHANICAL PEELING

Removes and cleans impurities via vacuum pressure (15 kPa to 80 kPa).

○ ELECTROPORATION + LED

Increases the penetration of ONA Moisturizing Serum through temporary pores in the cells created by high voltage pulses over short durations.

○ LOW FREQUENCY ELECTROTHERAPY

Modulated low-frequency (1kHz) contributes to muscle movement in a controlled manner, with the aim of toning the facial structure.




CLEAN!

TECHNOLOGICAL SYNERGY

Chemical peeling + mechanical

A special solution that is expelled with great power and then suctioned with our radial tip removes blackheads and impurities, all without damaging the skin.



VACUUM SUCTION



Vacuum cleaning and exfoliation



Deep cleaning




Effective penetration of solutions

Solutions: ON1, ON2 and ON3




ON1 AHA

Formula created from a blend of naturally derived alpha-hydroxy acids (AHA) of natural origin results in soothing and regenerating extracts, specially designed to exfoliate, clean, and deeply care for the skin.



ON2 BHA

Formula designed to purify and regenerate the skin from salicylic acid (BHA), willow extract, and other soothing and regenerating active ingredients, specially designed to treat clogged pores and skin imperfections. Ideal for the T-zone.



ON3 Anti-Aging Complex


Formula specially designed to restore the barrier function, deeply moisturize, and illuminate the skin, thanks to a blend of hyaluronic acid, regenerating peptides, and plant extracts with anti-aging effect.

Treatment Tips



The ON1, ON2 and ON3 solutions circulate radially through the single-use treatment tips, available in small and large sizes, to customize the treatment, depending on the type of skin or area to be treated.

This radial movement allows the pore to be deeply cleaned, without damaging the skin.



THREE
TWO
ON1

INFUSE AND TREAT!

TECHNOLOGICAL SYNERGY

Electroporation + LED

Increases and accelerates the penetration of the active ingredients of the ONA Moisturizing Serum up to 5mm deep, due to the temporary channels created by destabilizing the cell membrane, thanks to high-voltage pulses.

Through this process, the biological structure and cellular function are not damaged or altered.



Effective penetration of ONA moisturizing serum



Hydration of the skin

Effective phototherapy

Different wavelengths of LED lights create a synergy with electroporation, boosting results.



BLUE LED


Soothing effect on oily skin tendencies

RED LED

Stimulates the production of collagen, with great anti-aging effect

VIOLET LED

Combined effect (blue + red)



ONA Moisturizing Serum

Concentrated serum with a powerful anti-aging effect, thanks to its content of moisturizing and repairing plant extracts, medium molecular weight hyaluronic acid, and peptides that stimulate collagen production, which reinforce the skin's barrier function.

ONA Moisturizing Serum has been specifically formulated by INDIBA® with active ingredients of the highest quality. Its ingredients include:

- **Medium molecular weight hyaluronic acid**

Hyaluronic acid with a specific size that stimulates damaged and/or aged epithelial cells, providing smoother, more uniform, and softer skin.

- **Betaine**

With soothing properties, it strengthens the skin barrier by increasing cell integrity and cohesion, reducing redness, and protecting the skin from exposure to harmful substances.

- **Anti-aging peptides**

Small endogenous peptides derived from the proteolysis of the extracellular matrix, capable of regulating, repairing, and regenerating damaged and/or aged skin.



THREE
TWO
ONA

WORKOUT!

TOTAL ANTI-AGING EFFECT

Low-frequency electrotherapy + massage

Stimulates the facial musculature, contracting and relaxing the tissue repeatedly. It also reduces tension and improves nutrient delivery by increasing local blood circulation.



Low-frequency current (1 kHz)



Compact spheres with 360° rotation



Improves skin elasticity




Muscle massage by deep stimulation




Skin regeneration, muscle stimulation, and elimination of impurities



Promotes better circulation



MASSAGE



LOW FREQUENCY (1kHz)

Protective gel for electrotherapy ONA

Conductive gel specially formulated to maximize ONA results, safely and effectively.



THREE
TWO
ONA

THREE
TWO
ONA



Technical Specifications

Input power	220-240V, 50/60Hz
Input fuse	2 xT 10AL 250V
Power consumption	250 Watts
Weight	25 kg (including solutions)
Size	448mm x 438mm x 391mm (H x W X D)
Screen	7-inch LCD touchscreen
Cleansing Head	Level 1~7, Max 630 mmHg
Infuse Head	Level 0~10, Max 2.0 Vrms (@ 500 ohm)
Workout Head	Level 0~10, Max 3.8 Vrms (@ 500 ohm)

De Grift 20
7711 EJ, Nieuwleusen

www.laservision.nl | info@laservision.nl | 0529 - 428000 | [@laservision_instruments](https://www.instagram.com/laservision_instruments)